Indications for Use:

IneoGuard™ glove is a disposable sterile powder-free surgical glove made of synthetic elastomer that is intended to be worn by operating room personnel to protect a surgical wound from contamination. The glove contains radiation-attenuating particles that provided only a limited efficiency to reduce the exposure from harmful ionizing scattered or secondary radiation on the operator's hand during fluoroscopic procedures.

Caution:

These gloves are not to be used in or next to the primary X-Ray beam. The choice of a suitable glove must be made according to the activity being performed, the type of risk and the related environmental condition.

Radiation attenuation gloves offer limited protection to health care providers exposed to scattered radiation from patients during fluoroscopically guided procedures. The hands of the fluoroscope operator should only be placed in the Field of View when required by the goals of the procedure, whether a radiation attenuation glove is used or not.

This product is also classed as Personal Protective Equipment (PPE) according to European PPE Regulation 2016/425 and has been shown to comply with this Directive through the Harmonized European Standards EN 388:2016+A1:2018, EN1420:2020, EN421:2010 and EN374-1:2016, EN 16523-1:2015+A1:2018, EN374-2:2014, EN374-4:2013 and EN374-5:2016

PERFORMANCE AND LIMITATION OF USE – This product has been tested in accordance with EN 388:2016+A1:2018, EN1420:2020, EN421:2010, EN374-1:2016, EN 16523-1:2015+A1:2018, EN374-2:2014, EN374-4:2013, EN374-5:2016 and EN61331-1:2014 and achieved the following performance levels:

* EN388:2016+A1:2018: Levels of performance [abrasion 0 / cut 0 / tear 0 / puncture 0]
This glove does not protect against mechanical risks

* EN374-1:2016 and EN 16523-1:2015+A1:2018: Permeation levels are based on breakthrough times as follows:

Level 1 > 10mins	Level 3 > 60mins	Level 5 > 240mins
Level 2 > 30mins	Level 4 > 120mins	Level 6 > 480mins

EN ISO 374-1:2016 / Type B	Level	EN374-4:2013 % Degradation
40% Sodium hydroxide (K)	6	-24.6%
65% Nitric Acid (M)	3	8.1%
30% Hydrogen Peroxide (P)	5	-2.5%
37% Formaldehyde (T)	5	-14.1%

This information does not reflect the actual duration of protection in the workplace and the differentiation between mixtures and pure chemicals.

The chemical resistance has been assessed under laboratory conditions from samples of IneoGuard™ Model 1 taken from the palm only and relates only to the chemical tested. It can be different if the chemical is used in a mixture.

It is recommended to check that the gloves are suitable for the intended use because the conditions at the workplace may differ from the type test depending on temperature, abrasion and degradation.

When used, protective gloves may provide less resistance to the dangerous chemical due to changes in physical properties. Movements, snagging, rubbing, degradation caused by the chemical contact etc may reduce the actual use time significantly. For corrosive chemicals, degradation can be the most important factor to consider in selection of chemical resistant gloves.

Before usage, inspect the gloves for any defect or imperfections.
* EN374-2: Water-leak tested by Manufacturer
Level 2: AQL<1.5: PASS

* EN374-5: Resistance to penetration by blood-borne pathogens – test method using Phi-X174 bacteriophage: PASS

* EN421-2010 & EN61331-1 (2014) – Narrow Beam Configuration
Complies with 4.2.1; 4.2.2; 4.3 (Model 2 only); 4.4; 4.5 and 5.0

Lead Equivalency at 80kVp:

IneoGuard Model 1: 0.05mm
IneoGuard Model 2: 0.07mm

Average % attenuation values:

X-ray tube voltage (kV)	60	80	100	120
IneoGuard Model 1	52%	44%	40%	36%
IneoGuard Model 2	61%	54%	49%	45%

Narrow Beam Geometry. Sampling based on EN421, average on 4 locations and 2 gloves

Warning:

This glove does not protect against mechanical risks
These gloves contain elements of high atomic number and shall not be used to protect against radioelement of high beta energy so as to avoid bremsstrahlung

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