ineoGuard™

ENGLISH (EN)

ineoGuard™

INSTRUCTIONS FOR USE

User Notice
Warning: Read Carefully Before Use

This product complies with the European Medical Device Regulation 2017/745 (class IIa device) and is also certified as category III Personal Protective Equipment according to the European PPE Regulation 2016/425.

This product satisfies the applicable requirements of the following Harmonized European standards: EN 455-1:2020+A2:2024, EN 455-2:2024, EN 455-3:2023, EN 455-4:2009; EN 61331-1:2014; EN ISO 10993-5:2009+A11:2025; EN ISO 10993-10:2023 ; EN ISO 13485:2016+A11:2021, EN ISO 14971:2019+A11:2021, EN ISO 15223-1:2021+A1:2025, EN ISO 11607-2:2020+A1:2023, EN ISO 11137-1:2026, EN ISO 11137-2:2015+A1:2023, EN ISO 21420:2020+A1:2024, EN 421:2010 and EN ISO 374-1:2016+A1:2018, EN 16523-1:2015+A1:2018, EN ISO 374-2:2019, EN ISO 374-4:2019 and EN ISO 374-5:2016

Product Performance:

* **EN 455-1:** Freedom from hole, AQL=0.65
* **EN 455-2:** Comply with dimensions and physical properties requirements.
* **EN 455-3:** Absence of powder

* **EN ISO 374-1:2016+A1:2018 and EN 16523-1:2015+A1:2018:** Permeation levels are based on breakthrough times as follows:

Level 1 > 10mins	Level 3 > 60mins	Level 5 > 240mins
Level 2 > 30mins	Level 4 > 120mins	Level 6 > 480mins

EN ISO 374-1:2016+A1:2018 / Type B	Level	EN ISO 374-4:2019
40% Sodium hydroxide (K)	6	-36.6%
65% Nitric Acid (M)	3	4.7%
30% Hydrogen Peroxide (P)	5	-13.3%
37% Formaldehyde (T)	5	-43.8%

The chemical resistance has been assessed under laboratory conditions from samples of IneoGuard™ Model 1 taken from the palm only and relates only to the chemical tested; results may differ if the chemical is used in a mixture. EN ISO 374-4:2019 degradation levels indicate the change in puncture resistance after exposure to the challenge chemical.

The penetration resistance has been assessed under laboratory conditions and relates only to the tested specimen.

It is recommended to check that the gloves are suitable for the intended use because the conditions at the workplace may differ from the type test depending on temperature, abrasion and degradation.

When used, protective gloves may provide less resistance to the dangerous chemical due to changes in physical properties. Movements, snagging, rubbing, degradation caused by the chemical contact etc. may reduce the actual use time significantly. For corrosive chemicals, degradation can be the most important factor to consider in selection of chemical resistant gloves.

* **EN ISO 374-2: Air Leak (3.00kPa) & Water-leak tested. -** No Leaks Detected **: PASS**

* **EN ISO 374-5: Resistance to penetration by blood-borne pathogens –** test method using Phi-X174 bacteriophage**: PASS**

*** EN 421-2010 & EN 61331-1 (2014) – Narrow Beam Configuration**

Lead Equivalency at 80kVp:

IneoGuard Model 1: 0.05mm
IneoGuard Model 2: 0.07mm

Average % attenuation values:

X-ray tube voltage (kV)	60	80	100	120
IneoGuard Model 1	52%	44%	40%	36%
IneoGuard Model 2	61%	54%	49%	45%

Narrow Beam Geometry. Sampling based on EN421, average on 4 locations and 2 gloves. Typical performance level: all gloves with ≤3% uncertainty.

***ASTM 2547 & ASTM 7866**

X-ray tube voltage (kV)	60	80	100	120
IneoGuard Model 1	41%	34%	28%	24%
IneoGuard Model 2	50%	42%	36%	31%

Measured at the thinnest glove portion (cuff), average on 13 pieces/batch. Typical performance level: all gloves with ≤3% uncertainty.

Both models comply with min. attenuation requirements defined in ASTM 7866

Instructions for storage:

The glove must be stored in its original packaging at room temperature in a dry place away from sunlight.

Before use:

The gloves are available in sizes 5.5, 6.0, 6.5, 7.0, 7.5, 8.0, 8.5 and 9.0 and are supplied sterile (SAL=10^{-6} , ISO 11137).

ineoGuard™

INSTRUCTIONS FOR USE

User Notice
Warning: Read Carefully Before Use

Check the expiry date: The glove may only be used until the expiry date indicated on the packaging, provided that it has been stored properly in its intact packaging and under appropriate storage conditions until this date.

Prior to use, visually inspect the pouch to ensure its complete integrity. A pouch that is not intact (visible damage such as holes, cut, compromised sealing) must be discarded. Before usage, inspect the gloves for any defect or imperfections.

Donning and doffing instruction:

Don the glove carefully to avoid cuff tearing. For detailed donning and doffing instructions, please refer to: https://ineoguard.com/innovate-with-ineotech/. The used gloves shall be disposed of as clinical waste.

Recommended Use:

IneoGuard™ glove is a disposable sterile powder-free surgical glove made of synthetic elastomer that is intended to be worn by clinicians during fluoroscopic-guided procedures to protect a surgical wound from contamination. The glove contains radiation-attenuating particles that provided only a limited efficacy to reduce the exposure from harmful ionizing scattered or secondary radiation on the operator's hand during fluoroscopic procedures. Target groups are composed of healthcare providers performing fluoroscopic guided procedures as well as their patients. The glove is suitable for procedures performed under fluoroscopy where a surgical glove and limited level of attenuation to ionizing radiations are needed such as in: interventional radiology, interventional cardiology, neuroradiology, intensive care, endoscopy, urology, orthopedics and spine surgery.

After use:

Gloves are supplied as single use.

Caution:

The gloves do not protect against mechanical risks (abrasion, cuts, tears, piercing, etc..) and against heat and fire.

These gloves are not to be used in or next to the primary X-Ray beam. The choice of a suitable glove must be made according to the activity being performed, the type of risk and the related environmental condition.

Radiation attenuation gloves offer limited protection to clinicians exposed to scattered radiation from patients during fluoroscopically guided procedures. The hands of the clinician should only be placed in the Field of View when required by the goals of the procedure, whether a radiation attenuation glove is used or not.

These gloves contain elements of high atomic number and shall not be used to protect against radioelement of high beta energy so as to avoid bremsstrahlung.

Any serious incident that has occurred in relation to the device should be reported to IneoTech and the competent authority of the European Member State in which the user and/or patient is established.

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Euro Diffusion Medical SAS